

# Tris-EDTA 100X Solution

<b>Product</b>	Tris-EDTA 100X Solution
<b>Manufacturer</b>	Fisher Scientific
<b>Composition</b>	85% water, 12.1% tris (hydroxymethyl) aminomethane, 2.9% ethylenediamine tetraacetic acid, <1% hydrochloric acid; liquid
<b>Synonyms</b>	Tris(hydroxymethyl)aminomethane-EDTA solution
<b>Typical Uses</b>	Electrolyte solution component
<b>Hazards</b>	1 0 0 Not a hazardous substance or mixture.
<b>Storage</b>	Not stored at the nanoFAB.
<b>Disposal</b>	Dump in drain of Wet Deck 1A/1B and rinse copiously with DI water.
<b>Comments</b>	May be mixed with water and/or potassium chloride at Wet Deck 1A/1B.
<b>Status</b>	



## Check First

Please contact a nanoFAB staff member for uses not listed above.