

# Isopropyl alcohol (IPA)

Product	Isopropyl alcohol
Manufacturer	ThermoFisher Scientific
Composition	>95% 2-propanol
Synonyms	IPA, isopropanol, 2-propanol, propan-2-ol
Typical Uses	Organic solvent
Hazards	<div></div> <p>203</p> <p>Flammable liquids <b>CATEGORY 2</b>; Serious eye damage/eye irritation <b>CATEGORY 2</b>; Specific target organ toxicity – single exposure <b>CATEGORY 2</b> (kidney, liver); Specific target organ toxicity – single exposure <b>CATEGORY 3</b> (central nervous system)</p>
Storage	Cabinet 7 (flammables cabinet in Wet Aisle 2); underneath litho wet decks
Disposal	Solid: Organic waste / organic sharps waste containers in primary or secondary litho areas. Liquid: solvent waste bottles in litho wet decks or wet decks 1A/1B/2A. Cleanroom wipes lightly moistened with IPA (e.g., for cleaning items in gowning) may be disposed of in regular garbage bins.
Comments	Heavily used within the nanoFAB. IPA can be found in squeeze bottles with yellow lids throughout the centre.
Status	STOCKED



## Check First

Please contact a nanoFAB staff member for uses not listed above.

## Additional SDS Sheets

- [Fisher Scientific - Isopropanol](#)