




Acetone

Product	Acetone
Manufacturer	Fisher Scientific
Composition	100% acetone ; liquid
Synonyms	2-propanone, propanone, dimethyl ketone
Typical Uses	Organic solvent; cleaning; photoresist removal
Hazards	<div>    </div> <p>203</p> <p>Flammable liquids CATEGORY 2; Serious eye damage/eye irritation CATEGORY 2; Specific target organ toxicity – repeated exposure CATEGORY 2 (kidney, liver, spleen, blood); Specific target organ toxicity – single exposure CATEGORY 3 (central nervous system)</p>
Storage	Cabinet 7 (flammables cabinet in Wet Aisle 2); underneath litho wet decks
Disposal	Solid: Organic waste / organic sharps waste containers in primary or secondary litho areas. Liquid: solvent waste bottles in litho wet decks or wet decks 1A/1B/2A. Cleanroom wipes lightly moistened with acetone may be disposed of in regular garbage bins.
Comments	Heavily used within the nanoFAB. Acetone can be found in squeeze bottles with red lids throughout the centre.
Status	STOCKED



Check First

Please contact a nanoFAB staff member for uses not listed above.